

FUME HOODS

MOBILE DUCTLESS FUME HOOD

Filter & Low Air
Flow Alarm for Safety

Quick Filter
Change System


Chemical Resistant
And High Impact
Polycarbonate Sash

Highly
Chemical
Resistant
Worktop

Utility
Grommet

Ergonomic
Handtrap

High Quality
Lockable Castors



DLF-12M / DLF-15M

METHOD INNOVATIVE AND MODERN DUCTLESS FILTERING SYSTEM

SAFETY, CONVENIENCE AND AFFORDABILITY

Features

- 1) NFX-211 tested ductless filtration system from Europe.
- 2) Very low energy consumption, reducing operating costs.
- 3) Ready to use and no ducting need.
- 4) Mobile unit allows for flexibility and convenience.
- 5) Ergonomic handtraps for optimal user protection.
- 6) Large working space, almost full length of use.
- 7) Highly efficient and silent fan operation.
- 8) Audible and visual alarm to alert user of filter saturation and low air flow.
- 9) Convenient and easy to operate system.
- 10) Electrical grommets on both sides to allow for easy cable handling.
- 11) Filter view window to indicate filter model and correct placement of filter.



Utility Grommet



Ergonomic Handtrap



Filter & Air Flow Alarm



Easy Filter Replacement System

Technical Specification

Model	DLF-12M	DLF-15M
Size (mm) w/o base:	1220(L) x 1180(H) x 800(D)	1500(L) x 1180(H) x 800(D)
Structure:	Steel construction epoxy baked finish with three sided tempered glass & flip open polycarbonate front sash.	
Internal work size:	1148(L) x 950(H) x 773(D)	1428(L) x 950(H) x 773(D)
Size (mm) with base:	1220(L) x 2040(H) x 800(D)	1500(L) x 2040(H) x 800(D)
Workshop:	Epoxy Resin Worktop / Compact Phenolic Laminate.	
Head System:	DLF-TH301SFA (Single Filter)	DLF-TH501SFA (Dual Filters)
Filters:	Active Coal Filters for Chemicals Fumes	
Voltage:	220-240V	

*Please ask us for other filter types



METHOD ENTERPRISE SDN. BHD. (195209-M)

Tel : +6 03 5122 1818 Fax : +6 03 5122 8899

Email : sales@method.com.my / method_enterprise@yahoo.com

Website : www.method.com.my / www.lab-manufacturer.com